

# Automatic Dose Optimization for Radiotherapy

Paul Dubois

TheraPanacea  
MICS, CentraleSupélec  
Institut du Cancer de Montpellier


# Re-inventing cancer treatment with Artificial Intelligence



**THERAPANACEA**


Reinventing cancer care through AI

# Annotation



**THERAPANACEA**  
Reinventing cancer care through AI

# Annotation



**THERAPANACEA**  
Reinventing cancer care through AI



# Annotate

## AI-powered contouring tool



In clinical use since 2019



150 + structures (OAR and LN) on CT

→ Sub structures of the heart

→ SBRT Thorax

70 + structures on MR



3-4 minutes contouring



Plug and play



**THERAPANACEA**  
Reinventing cancer care through AI

# MR Box

## Synthetic CT for MR-guided workflow



Automatic delineation on MR



Generate a pseudo-CT from MR images

- Avoid registrations errors
- Accelerate adaptive routines
- Less machine time



Brain T1, Pelvis T2, Pelvis/Abdo Truefisp (ViewRay)  
And more under development



**THERAPANACEA**  
Reinventing cancer care through AI

ART-Plan™ AdaptBox

# Advantages

Accelerate, automatize and simplify adaptive radiation therapy.



Reduce the burdens of multiple manual iterations.



Save time and optimize the resources of your department.



Take earlier and easier re-planning decision.

# Cancer treatments

Surgery



Chemotherapy



- + : Safe (little damage to healthy tissues)
- : Tumor needs to be localized & accessible



- : Heavy medicine on all the body
- + : Tumor does **not** need to be localized




# Cancer treatments

Surgery



Radiotherapy



Chemotherapy




+ : Safe  
- : Tumor needs to be localized


+ : Relatively safe (most tissues are spared)  
- : Tumor needs to be (relatively) localized

Chemotherapy on all the body  
Tumor does **not** need to be localized


# Multi-Leaf Collimator



# Spreading rays




# Dose Spread







# Irradiation strategy



# Delivering irradiation



Ideal heatmap




Discretized heatmap

# Delivering discretized dose



*dose to deliver*



# Problem statement



% of the volume	Name	Dose
<b>PTV constraints:</b>		
$\leq 1\%$	of Ext_Tumor receives	$\geq 70\text{Gy}$
$\geq 95\%$	of Tumor receives	$\geq 63\text{Gy}$
<b>OAR constraints:</b>		
$\leq 1\%$	of Spinal_Cord receives	$\geq 43\text{Gy}$
$\leq 15\%$	of Heart receives	$\geq 30\text{Gy}$
$\leq 20\%$	of Esophagus receives	$\geq 10\text{Gy}$
$\leq 2\%$	of Lt_Lung receives	$\geq 20\text{Gy}$
$\leq 8\%$	of Lt_Lung receives	$\geq 10\text{Gy}$

List of “objectives”



# Mathematical formulation

Parameters to optimize

$$x_{i,j}^{(\theta)} \text{ with } 1 \leq i, j \leq N_{leaf\_pairs}$$
$$\theta \in \Theta \text{ (the list of angles chosen for this case)}$$

Parameters constraints

We cannot send negative dose, i.e.:

$$x_{i,j}^{(\theta)} \geq 0$$

however, the angles are completely free.

Dose calculation

The dose on each voxel of the body  $y_v$  is approximated by:

$$y = \sum_{\theta \in \Theta} L^{(\theta)} x^{(\theta)}$$

We want to find the best dose distribution  $y$ .

# Mathematical formulation

Ideal problem

$$\begin{aligned} & \text{minimize} \sum_{s \in \mathcal{O}} \frac{1}{|\mathcal{V}_s|} \sum_{v \in \mathcal{V}_s} y_v^2 \\ & \text{s.t. DVH are respected, } y = Lx \text{ and } x \geq 0 \end{aligned}$$

% of the volume	Name	Dose
<b>PTV constraints:</b>		
≤ 1%	of Ext_Tumor receives	≥ 70Gy
≥ 95%	of Tumor receives	≥ 63Gy
<b>OAR constraints:</b>		
≤ 1%	of Spinal_Cord receives	≥ 43Gy
≤ 15%	of Heart receives	≥ 30Gy
≤ 20%	of Esophagus receives	≥ 10Gy
≤ 2%	of Lt_Lung receives	≥ 20Gy
≤ 8%	of Lt_Lung receives	≥ 10Gy

Objective function



$$f_s(y) = \sum_{v \in \mathcal{V}_s} \frac{1}{|\mathcal{V}_s|} \left( \overline{w}_s (y_v - \overline{t}_s)_+^p + \underline{w}_s (-y_v + \underline{t}_s)_+^p \right)$$


$$f(y) = \sum_{s \in \mathcal{S}} f_s(y)$$

*if constraint is not satisfied for structure  $s$ , else  $f_s(y) = 0$*

# Current Workflow

## Patient data

- CT scan
- OARs & PTVs contours
- doctors' objectives



## Treatment plan

- Fluences
- Leaf movements
- Dose per organ



## Input from dosimetrist

- Weights of each objective


# Automatization (*work in progress*)

## Patient data


- CT scan
- OARs & PTVs contours
- doctors' objectives

## Treatment plan

- Fluences
- Leaf movements
- Dose per organ



Learning based on previous cases







# Thank You

